

Product Datasheet

F box and leucine-rich-repeat gene 4 Antibody (orb43386)

Description	F box and leucine-rich-repeat gene 4 Antibody
Species/Host	Rabbit
Reactivity	Drosophila, Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9VY46
Expiration Date	12 months from date of receipt.