

Product Datasheet

CD79a antibody (RPE) (orb433867)

Description	Mouse monoclonal antibody against CD79a conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Synthetic peptide corresponding to 202-216 amino acid sequence of human mb-1
Target	CD79a
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.