

Product Datasheet

CD79a antibody (orb433869)

Description	Mouse monoclonal antibody to CD79a
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P, WB
Immunogen	Synthetic peptide corresponding to 202-216 amino acid sequence of human mb-1
Target	CD79a
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.