

Product Datasheet

PLVAP antibody (Biotin) (orb433871)

Description	Rat monoclonal antibody against PLVAP conjugated to Biotin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	IH-Fr
Target	PLVAP
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.