

## Product Datasheet

# FELINE IMMUNODEFICIENCY VIRUS gp95 antibody (orb433882)

|                                |                                                                       |
|--------------------------------|-----------------------------------------------------------------------|
| <b>Description</b>             | Mouse monoclonal antibody to<br>FELINE IMMUNODEFICIENCY VIRUS<br>gp95 |
| <b>Species/Host</b>            | Mouse                                                                 |
| <b>Reactivity</b>              | Virus                                                                 |
| <b>Conjugation</b>             | Unconjugated                                                          |
| <b>Tested<br/>Applications</b> | FA, FACS, WB                                                          |
| <b>Target</b>                  | FELINE IMMUNODEFICIENCY VIRUS<br>gp95                                 |
| <b>Preservatives</b>           | 0.1% Sodium Azide (NaN <sub>3</sub> ).<br>Phosphate buffered saline   |
| <b>Form/Appearance</b>         | Purified IgG - liquid                                                 |
| <b>Concentration</b>           | 1.0mg/ml                                                              |
| <b>Storage</b>                 | +4°C, -20°C if preferred                                              |
| <b>Note</b>                    | For research use only                                                 |
| <b>Application notes</b>       | V3 REGION                                                             |
| <b>Isotype</b>                 | IgG1                                                                  |
| <b>Clonality</b>               | Monoclonal                                                            |
| <b>Purity</b>                  | Purified                                                              |
| <b>Expiration Date</b>         | 12 months from date of receipt.                                       |