

Product Datasheet

CD107b antibody (orb433887)

Description	Mouse monoclonal antibody to CD107b
Species/Host	Mouse
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IP, WB
Immunogen	MDCK (Madin-Darby Canine Kidney) cells
Target	CD107b
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.