

Product Datasheet

ASPERGILLUS antibody (orb433901)

Description	Mouse monoclonal antibody to ASPERGILLUS
Species/Host	Mouse
Reactivity	Fungi
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, IP, WB
Immunogen	Wall fraction (WF) of Aspergillus fumigatus
Target	ASPERGILLUS
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgM - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	SPP
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.