

Product Datasheet

CD172g antibody (RPE) (orb433911)

Description	Mouse monoclonal antibody against CD172g conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Purified recombinant SIRPgamma.CD4.
Target	CD172g
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.