

Product Datasheet

SUPEROXIDE DISMUTASE antibody (orb433927)

Description	Mouse monoclonal antibody to SUPEROXIDE DISMUTASE
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr, IHC-P
Immunogen	Superoxide dismutase enzyme prepared from bovine erythrocytes.
Target	SUPEROXIDE DISMUTASE
Preservatives	0.1% Sodium Azide (NaN ₃) 0.1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(Cu-Zn)
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.