

Product Datasheet

CD117 antibody (orb433930)

Description	Mouse monoclonal antibody to CD117
Species/Host	Mouse
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP, WB
Immunogen	Porcine bone marrow cells (BMC).
Target	CD117
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.