

Product Datasheet

C1q antibody (orb433935)

Description	Mouse monoclonal antibody to C1q
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF, IHC-Fr, WB
Immunogen	Globular head domain of C1q, purified from human plasma.
Target	C1q
Preservatives	0.09% Sodium Azide (NaN ₃). Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.