

Product Datasheet

C4c antibody (orb433940)

Description	Mouse monoclonal antibody to C4c
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr, WB
Immunogen	Purified human C4c.
Target	C4c
Preservatives	0.09% Sodium Azide (NaN ₃). Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.