

Product Datasheet

C9 antibody (orb433993)

Description	Mouse monoclonal antibody to C9
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified human C9.
Target	C9
Preservatives	0.1% Sodium Azide (NaN ₃). Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.