

Product Datasheet

COMPLEMENT FACTOR B antibody (orb433994)

Description	Mouse monoclonal antibody to COMPLEMENT FACTOR B
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, IHC-Fr, WB
Immunogen	Purified human complement factor b.
Target	COMPLEMENT FACTOR B
Preservatives	0.1% Sodium Azide (NaN ₃). Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(Ba FRAGMENT)
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.