

Product Datasheet

C3d antibody (orb433995)

Description	Mouse monoclonal antibody to C3d
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr
Immunogen	Purified human C3d (neoantigen).
Target	C3d
Preservatives	0.1% Sodium Azide (NaN ₃). Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.