

Product Datasheet

CD137 antibody (orb433999)

Description	Rat monoclonal antibody to CD137
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Recombinant human CD137-Fc.
Target	CD137
Preservatives	0.1% Sodium Azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.