

Product Datasheet

EMR1 antibody (RPE) (orb434021)

Description	Hamster monoclonal antibody against EMR1 conjugated to RPE
Species/Host	Hamster
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Fibroblast cell line ARH102 expressing CD97-EMR1 (EGF1-6)
Target	EMR1
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.