

Product Datasheet

CD257 antibody (orb434024)

Description	Mouse monoclonal antibody to CD257
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant 17 kDa polypeptide corresponding to the TNF-like portion of the extracellular domain of human CD257, expressed in E. coli.
Target	CD257
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0mg/ml
Storage	-20°C only
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.