

Product Datasheet

CD3 antibody (orb434048)

Description	Hamster monoclonal antibody to CD3
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IP, WB
Immunogen	H-2Kb - specific mouse cytotoxic T lymphocyte clone BM10-37.
Target	CD3
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.