

Product Datasheet

HEXOKINASE 1 antibody (orb434067)

Description	Mouse monoclonal antibody to HEXOKINASE 1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant human hexokinase 1 (amino acids 1-917) purified from Escherichia coli.
Target	HEXOKINASE 1
Preservatives	0.02% Sodium Azide (NaN ₃) 10% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.