

Product Datasheet

IRF5 antibody (orb434070)

Description	Mouse monoclonal antibody to IRF5
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, IP, WB
Immunogen	Histidine-tagged recombinant peptide corresponding to amino acids 176-240 of human IRF5, purified from E. coli.
Target	IRF5
Preservatives	0.02% Sodium Azide (NaN ₃) 10% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.