

Product Datasheet

ASTROVIRUS antibody (orb434083)

Description	Mouse monoclonal antibody to ASTROVIRUS
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IP
Immunogen	Astrovirus type 2 expressed in LLC MK2 cells.
Target	ASTROVIRUS
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.