

Product Datasheet

CHLAMYDIA LPS antibody (orb434085)

Description	Mouse monoclonal antibody to CHLAMYDIA LPS
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Elementary bodies from C. trachomatis strain SAF2.
Target	CHLAMYDIA LPS
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.