

Product Datasheet

CD44 antibody (orb434092)

Description	Mouse monoclonal antibody to CD44
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr, IHC-P, IP, WB
Immunogen	Stimulated human leucocytes.
Target	CD44
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.