

Product Datasheet

CD172a antibody (Biotin) (orb434126)

Description	Mouse monoclonal antibody against CD172a conjugated to Biotin
Species/Host	Mouse
Reactivity	Rat
Conjugation	Biotin
Tested Applications	FC, IH-Fr, IHC-P
Immunogen	Rat peritoneal macrophages.
Target	CD172a
Preservatives	0.09% Sodium Azide 1.0% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.