

Product Datasheet

CD172a antibody (orb434130)

Description	Mouse monoclonal antibody to CD172a
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P, WB
Immunogen	Rat peritoneal macrophages.
Target	CD172a
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.