

## Product Datasheet

### FIBRINOPEPTIDE A antibody (orb434138)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to FIBRINOPEPTIDE A                    |
| <b>Species/Host</b>        | Mouse                                                            |
| <b>Reactivity</b>          | Human                                                            |
| <b>Conjugation</b>         | Unconjugated                                                     |
| <b>Tested Applications</b> | ELISA                                                            |
| <b>Immunogen</b>           | Synthetic fibrinopeptide A conjugated to a carrier protein.      |
| <b>Target</b>              | FIBRINOPEPTIDE A                                                 |
| <b>Preservatives</b>       | 0.1% Sodium Azide (NaN <sub>3</sub> ). Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                            |
| <b>Concentration</b>       | 1.0mg/ml                                                         |
| <b>Storage</b>             | +4°C, -20°C if preferred                                         |
| <b>Note</b>                | For research use only                                            |
| <b>Isotype</b>             | IgG2a                                                            |
| <b>Clonality</b>           | Monoclonal                                                       |
| <b>Purity</b>              | Purified                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |