

## Product Datasheet

### IgG2a antibody (orb434166)

|                            |                                              |
|----------------------------|----------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to IgG2a           |
| <b>Species/Host</b>        | Mouse                                        |
| <b>Reactivity</b>          | Rat                                          |
| <b>Conjugation</b>         | Unconjugated                                 |
| <b>Tested Applications</b> | ELISA                                        |
| <b>Immunogen</b>           | purified rat IgG2a.                          |
| <b>Target</b>              | IgG2a                                        |
| <b>Preservatives</b>       | 0.1% Sodium Azide. Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                        |
| <b>Concentration</b>       | 1 mg/ml                                      |
| <b>Storage</b>             | +4°C, -20°C if preferred                     |
| <b>Note</b>                | For research use only                        |
| <b>Application notes</b>   | HEAVY CHAIN                                  |
| <b>Isotype</b>             | IgG1                                         |
| <b>Clonality</b>           | Monoclonal                                   |
| <b>Purity</b>              | Purified                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.              |