

Product Datasheet

CD11c antibody (orb434174)

Description	Mouse monoclonal antibody to CD11c
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FACS, IHC-Fr
Immunogen	Recombinant human CD11c.
Target	CD11c
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.