

Product Datasheet

PAPILLOMAVIRUS 16 ONCOPROTEIN E7 antibody (orb434198)

Description	Mouse monoclonal antibody to PAPILLOMAVIRUS 16 ONCOPROTEIN E7
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Oncoprotein E7 from human papillomavirus 16.
Target	PAPILLOMAVIRUS 16 ONCOPROTEIN E7
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.