

Product Datasheet

SALMONELLA PARATYPHI A LPS antibody (orb434208)

Description	Mouse monoclonal antibody to SALMONELLA PARATYPHI A LPS
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	DOT, ELISA, ICC, IF
Immunogen	Heat inactivated cells of Salmonella paratyphi A.
Target	SALMONELLA PARATYPHI A LPS
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	GROUP A 0-2 ANTIGEN
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.