

Product Datasheet

CD86 antibody (RPE) (orb434250)

Description	Mouse monoclonal antibody against CD86 conjugated to RPE
Species/Host	Mouse
Reactivity	Rat
Conjugation	RPE
Tested Applications	FC
Immunogen	HTLV-1 transformed Lewis-S1 cells.
Target	CD86
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.