

Product Datasheet

CA125 antibody (orb434254)

Description	Mouse monoclonal antibody to CA125
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	CA125 antigen purified from human ovarian carcinoma.
Target	CA125
Preservatives	0.1% Sodium Azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.