

Product Datasheet

PROSTATE SPECIFIC MEMBRANE ANTIGEN antibody (orb434263)

Description	Mouse monoclonal antibody to PROSTATE SPECIFIC MEMBRANE ANTIGEN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant fragment corresponding to amino acids 44-750 of human GCPIL.
Target	PROSTATE SPECIFIC MEMBRANE ANTIGEN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.