

Product Datasheet

V5-TAG antibody (FITC) (orb434274)

Description	Mouse monoclonal antibody against V5-TAG conjugated to FITC
Species/Host	Mouse
Reactivity	Virus
Conjugation	FITC
Tested Applications	FACS, ICC, IF
Immunogen	A recombinant peptide corresponding to the PK tag fused to a class II - restricted T helper cell epitope.
Target	V5-TAG
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.