

Product Datasheet

CHIMERIC HUMAN IgE ANTI NP antibody (orb434287)

Description	Mouse monoclonal antibody to CHIMERIC HUMAN IgE ANTI NP
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA
Immunogen	Hapten, 4-hydroxy-3-nitrophenylacetyl (NP).
Target	CHIMERIC HUMAN IgE ANTI NP
Preservatives	0.09% Sodium Azide. None present
Form/Appearance	Tissue Culture Supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgE
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.