

Product Datasheet

FSH antibody (orb434290)

Description	Mouse monoclonal antibody to FSH
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr, IHC-P
Immunogen	Human FSH.
Target	FSH
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	BETA 2
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.