

Product Datasheet

CD45RA antibody (orb434292)

Description	Mouse monoclonal antibody to CD45RA
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P
Immunogen	Purified Rat spleen L-CA
Target	CD45RA
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(B CELLS ONLY)
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.