

Product Datasheet

CD169 antibody (orb434308)

Description	Mouse monoclonal antibody to CD169
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IP
Immunogen	Rat Spleen cell homogenate
Target	CD169
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.