

Product Datasheet

CD8 antibody (orb434315)

Description	Rat monoclonal antibody to CD8
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr
Immunogen	Mouse L cells transfected with human CD8 gene.
Target	CD8
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.