

Product Datasheet

ACTIN F TYPE antibody (orb434318)

Description	Mouse monoclonal antibody to ACTIN F TYPE
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr
Immunogen	Human monocytes and U937 cell line.
Target	ACTIN F TYPE
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgM - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.