

Product Datasheet

HERPES SIMPLEX VIRUS 2 gG antibody (orb434326)

Description	Mouse monoclonal antibody to HERPES SIMPLEX VIRUS 2 gG
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IP, WB
Immunogen	HSV-2 strain 25766.
Target	HERPES SIMPLEX VIRUS 2 gG
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.