

Product Datasheet

CD45 antibody (RPE) (orb434333)

Description	Mouse monoclonal antibody against CD45 conjugated to RPE
Species/Host	Mouse
Reactivity	Rat
Conjugation	RPE
Tested Applications	FC
Immunogen	Rat thymocyte membrane glycoproteins.
Target	CD45
Preservatives	0.09% Sodium Azide1%Bovine Serum Albumin5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.