

Product Datasheet

TCR ALPHA/BETA antibody (FITC) (orb434340)

Description	Mouse monoclonal antibody against TCR ALPHA/BETA conjugated to FITC
Species/Host	Mouse
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	Rat T Blasts and Erythrocytes.
Target	TCR ALPHA/BETA
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.