

Product Datasheet

CD106 antibody (orb434366)

Description	Mouse monoclonal antibody to CD106
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr
Immunogen	Rat CD106-transfected L5178Y cells.
Target	CD106
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.