

Product Datasheet

CD3 antibody (orb434391)

Description	Mouse monoclonal antibody to CD3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr
Immunogen	Human infant thymocytes and lymphocytes from a patient with Sezary Syndrome.
Target	CD3
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.