

Product Datasheet

CD3 antibody (RPE) (orb434393)

Description	Mouse monoclonal antibody against CD3 conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Human infant thymocytes and lymphocytes from a patient with Sezary Syndrome.
Target	CD3
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.