

Product Datasheet

DECTIN-1 antibody (RPE) (orb434419)

Description	Mouse monoclonal antibody against DECTIN-1 conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	NIH3T3 cells expressing full length Beta-glucan receptor.
Target	DECTIN-1
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.