

Product Datasheet

S-TAG antibody (orb434424)

Description	Mouse monoclonal antibody to S-TAG
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, Purification
Immunogen	S-tag recombinant protein.
Target	S-TAG
Preservatives	None present. Borate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.