

## Product Datasheet

### FACTOR VIII antibody (orb434435)

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to FACTOR VIII                          |
| <b>Species/Host</b>        | Mouse                                                             |
| <b>Reactivity</b>          | Human                                                             |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Tested Applications</b> | ELISA, RIA, WB                                                    |
| <b>Immunogen</b>           | Affinity purified human Factor VIII.                              |
| <b>Target</b>              | FACTOR VIII                                                       |
| <b>Preservatives</b>       | 0.09% Sodium Azide (NaN <sub>3</sub> ). Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                             |
| <b>Concentration</b>       | 1.0mg/ml                                                          |
| <b>Storage</b>             | +4°C, -20°C if preferred                                          |
| <b>Note</b>                | For research use only                                             |
| <b>Isotype</b>             | IgG1                                                              |
| <b>Clonality</b>           | Monoclonal                                                        |
| <b>Purity</b>              | Purified                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |